

Quality Control for Genomic DNA Samples

Purpose: Using Qsep100 to evaluate the quality of genomic DNA samples.

Introduction: More than 1000 genomic DNA samples were checked for their integrity status in order to evaluate their quality for downstream qPCR analysis..

Advantages: For the quality control of genomic DNA samples, Qsep100 provides a rapid and convenient examine system with excellent reproducibility over the course of 1000 runs and ability to distinguish between the intact and degraded samples. Qsep100 provides crucial information, which samples can be used for downstream qPCR processing.

